



**New level of comfort and safety**

Affordable price, very early (VEM) mobilization, antibacterial technology, easier disinfection, alarm systems

Electric hospital bed with VEM functions

**Famed NEXO**



Electric hospital bed, ICU bed

# Famed NEXO

**THE COMFORT OF AN ELECTRIC BED  
AFFORDABLE PRICE  
MODERN APPROACH TO REHABILITATION  
OPTIMIZED FOR A LIMITED NUMBER OF HOSPITAL STAFF**



Famed NEXO is an electric hospital bed that provides comfort during hospital stay, both in general departments and in intensive care (ICU). The multitude of functionalities, technological solutions, high durability, easier disinfection and built-in rehabilitation guidelines have been combined with an affordable price.



## COMFORT

- Easy standing up and sitting on the bed, full safety - minimum height 390 mm.
- A multitude of available control options and alarm handling panel.
- Pre-programmed most commonly used positions, incl. examination position, vascular positions and disinfection positions for decontamination tunnels.
- Recovery support - indicators informing about the optimal inclination, VEM, boards that can be used during rehabilitation.
- Easy to disinfect, gently lowerable, antibacterial **SoftDropPlus™** side rails that can be lowered below the bed, operated with one hand.
- **USB sockets (2.1A)** for charging mobile devices makes the stay in the hospital more pleasant.



## SAFETY

- Helping hand in the fight against infections - plastic elements and varnish were made in the bacteriostatic **pSilver™** technology.
- High side rails along the entire length of the bed, **fully compliant with the medical requirements** of EN 60601-2-52.
- Quick removal of the head board with the possibility of using it as an **additional CPR support**.
- **Extended steel reinforcements**, a lockable board and additional stainless steel stabilizers allow for easy rolling of the bed and safe use of reverse-Trendelenburg function.
- Easily accessible CPR levers - even when the side rails are lowered, one on each side of the bed.
- The buttons activating the controls in the side rails and the panel prevent accidental change of bed position.
- Unlocked wheels alarm and patient leaving the bed alarm.



## DISINFECTION

- Fewer elements of the bed pallet, reduced spaces between panels, **chassis and moving parts covered in 91% allow to shorten the disinfection time**.
- Plastic and metal elements (varnish) were made with the addition of silver nanoparticles in the **pSilver™** technology.
- To improve the washing process, choose an IPX6 bed - protected against the ingress of liquids, also in beds decontamination stations.
- NEXO has a dedicated position for washing in BDS - a guarantee of proper drying of the bed and reduction of handling time during washing.
- The control cables are hidden in the side rails mechanisms, which reduces the number of hard-to-reach places during manual washing.



## ECONOMY

- **High versatility** - long-term, palliative and bariatric care units, the possibility of selecting options for ICU care.
- **pSilver™** bacteriostatic protection reduces the risk of HAI (Hospital Acquired Infections).
- Durable bed pallet filling (HPL) reduces the risk of losing or destroying individual elements and also **shortens maintenance and cleaning time**.
- **IPX6 version available** with solutions dedicated to decontamination tunnels.
- Moving parts and electronic elements **easily accessible to service teams**, the lack of non-standard tools required for repairs facilitates the work of internal services (hospital technical teams).



## NEXO FACILITATES THE WORK OF THE HOSPITAL STAFF

- Central wheel lock - keep patient safe with minimum effort.
- Foot controller - saves time and allows you to combine various activities.
- Function locks - they guarantee that the patient will not change the bed settings against the recommendations.
- Nurse Call system - button in bed (buzzer) or the possibility of connecting to the central hospital system (optional).
- One touch positions - reduce complicated combinations with buttons that set the bed to frequently used, useful positions.
- Easy disinfection - choose the variant of the bed for BDS or use the reduced time of manual disinfection thanks to the covering of the moving parts and the pSilver™ anti-bacterial silver technology.
- VEM (Very Early Mobilization) functions - very early mobilization of the patient, also in the ICU. Faster recovery of the patient.
- Side rails along the entire length of the bed - ensure the patient's safety regardless of the position of the NEXO. SoftDropPlus™ side rails can be lowered below the mattress platform level, very close to the structure of the bed - plenty of space around the bed, faster recovery of the patient without the risk of trapping the body or breaking the side rails
- The 5th wheel (option) and the rigid boards with extended reinforcements make it easier to roll the bed and overcome inclines or sharp turns.



# COMFORT

- A** **VEM central panel.** Preset positions (examination, emergency, Fowler, cardiology chair, CPR, two vascular positions), three VEM (Very Early Mobilization) settings, activation button, flashlight, locks, mains disconnection indicator, battery indicator.
- B** **Control for patient.** Activation button, position adjustments, under bed light, access to the nurse call system. Adapted to the visually impaired (contrasting icons and protrusions describing specific functions).
- C** **Control for personnel.** Activation button, controlling positions of the bed (including longitudinal tilt), highlighted emergency position.
- D** **VEM mobilization buttons.** They enable very early mobilization of the patient - helps to get out of bed, the patient rests on the railing, increasing the height with the button.
- E** **Foot control.** It allows you to adjust the height and at the same time work with the patient and assist in changing the position in the bed.



**Preprogrammed positions.** You can easily set useful and frequently used bed positions with a single button and save time.



**As many as 9 grip points** in the side rails, it allows you to conveniently and safely change positions and get out of bed.



**USB sockets.** Possibility to charge personal devices. Equipped with a sealed cap. Maximum current up to 2.1A.



**Optimal angle indicator.** Indicators help to position the patient in the prevention of ventilator pneumonia (VAP).



**Double autoregression.** Reduction of excessive pressure in a sitting position in the area of the lumbar spine and shin.



**Low position.** It allows for safe sleep and comfortable sitting on the bed and getting up. It also allows for easy transfer to a wheelchair.



## ONLY 35% HOSPITALS MEET THE STANDARDS OF COMFORT AND SAFETY<sup>1</sup>

We have created the Famed NEXO bed to enable hospitals to improve those standards and quickly raise these statistics



Although the statistics are improving (in 2015 the figure was 26%), still as many as 35% of examined hospitals do not meet the standards defined by an independent organization supervising the health and safety of patients, Leapfrog.

### The standards concern primarily:

- percentage of patients who get bedsores in the hospital,
- percentage of patients who had an incident of falling out of bed or other medical incident.



# SAFETY



**SoftDropPlus™ side rails along the entire length of the bed.** Gas springs ensure a quiet and smooth lowering motion. The constant distance between the elements ensures comfortable and safe usage, preventing accidental trapping of hands. The movement can be initiated with one hand, while the other hand can be used to protect the patient.



**Use the head board as a CPR support.** The unique design of the SoftDropPlus™ side rails gives an unprecedented amount of space for a resuscitation action or a mattress. Additionally, the quickly removable, stable head board can be used as additional support under the patient in the case of CPR action.



**CPR levers are easily accessible on both sides of the bed.** Convenient to use, easy to spot levers are available on both sides of the bed, regardless of the position of the side rails - in full range of motion of the backrest. The proximity to the patient makes it possible to check his condition on an ongoing basis when lowering the backrest in a life-threatening situation.



**use the best reverse-Trendelenburg angle on the market.** Foot board with lock, 18.5° tilt and extended, large diameter reinforcements allow you to safely move the bed and perform exercises and elements of upright standing therapy. Tilts also help to move the patient on the mattress.



**Controls with activation systems.** They reduce the risk of accidental change of position by the patient or the staff (during washing or handling). Only pressing the green activation button enables the work of the controllers.



**Central wheel lock.** The levers on both sides of the bed allow easy and quick blocking of the bed - they can also be operated from the front. Clear color coding and acoustic confirmation of the lever position make it easier to recognize the brake position.

# NEXO FOR DECONTAMINATION STATIONS

The bed variant in the IPX6 standard enriched with solutions that make the bed resistant to the conditions of BDS.

Famed Żywiec beds have been used in decontamination stations for years. Now after gathering experience with washing in the conditions of BDS we go one step further by introducing chemical protection. The NEXO bed for decontamination tunnels was made with attention to detail appropriate to the conditions that go with this disinfection method.

The entire structure, including the inside of the steel profiles, is protected using the modern CDC method (Cathodic dip coating). Protecting the bed after it has been subjected to metal processing guarantees high strength of welds, connections or holes intended for the outflow of cleaning fluids. Dedicated washing position, controllers prepared for washing in BDS, modifications and stainless details allow for trouble-free operation.



## UNIQUE HELPING HAND IN THE PREVENTION OF HAI

Study shows that hospital acquired infections (HAI) are a major threat for both health care providers and the patients<sup>2</sup>. Increased length of stay, rising costs of patient treatment, risk of lawsuits are a big danger for financial situation of the hospital. This is also an inconvenience and an increased risk of hospitalization for the patient. WHO guidelines on hand hygiene in health care clearly states that hands are the most common route of infection<sup>3</sup> - proper hygiene can be supported by technology to limit the bacteria spread. Equipping your hospital with the Famed NEXO bed with pSilver™ technology is one of the best ways to support the HAI prevention program in your facility. The use of pSilver™ technology in places most often touched by the patient and hospital staff reduces the multiplication of dangerous bacteria, fungi or viruses. All plastic and painted elements are made with custom compound with silver nanoparticles.





# DISINFECTION



**Bacteriostatic pSilver™ technology**  
Plastic elements contain silver nanoparticles - this is a lifetime protection and a guarantee of less frequent reproduction of viruses and bacteria. We use a similar antibacterial additive in the varnish powder. Famed Żywiec has been offering and refining this technology for several years.



**Smart-thinking: panels filling the mattress platform.** The reduced number of elements that staff have to disassemble for cleaning and the reduced number and size of spaces between the HPL elements means reduced cleaning and disinfection times. The mattress platform cover as much as 91% of the surface, protecting the actuators and pantograph chassis against contamination.



**Optimal cables routing.** A bed of this class has not only power cables, but also many control and grounding cables. The main one, connecting the bed with side rails has been hidden in the side rails mechanisms. The cables from the control panel are also covered, thanks to which they do not get dirty and are easily accessible for service teams.



**Dedicated version for automatic washing.** The version for washing stations is not only the IPX6 standard, it is also a completely new technology of protection and hundreds of hours of testing in laboratories and washing chambers.



**Dedicated button - washing position.** At the push of a button, all sections of the bed and the entire bed are set in a dedicated position for washing in BDS. Thanks to this, the bed is not exposed to damage, and the drying process runs smoothly.



**Tested resistance to disinfecting agents.** The more and more aggressive fight against germs means the constant development of disinfectants. Famed Żywiec, as one of the few manufacturers of medical equipment, has its own laboratory that tests every available agent. We deliver our products with a compatibility list and cleaning recommendations.

Antibacterial pSilver™ technology

Head board acting as a support for CPR

HPL transparent to X-rays (MAE 0.39 mm Al)

Bright, contrasting colors on the panel and protrusions on the icons dedicated for visually impaired.

Side rails along the entire length of the bed

Control panels dedicated to patients and staff

Foot board lock to prevent it from falling out during transport and enabling VEM therapy and rehabilitation

USB sockets (2.1A) - option

Possibility to personalize the bed

Lowest height of 390 mm

Plastic bumpers (less damage for medical equipment and walls)

Hooks with position adjustment

Easy access to CPR levers

Innovative VEM panel, a multitude of quickly available positions

Anatomical hand grip and very early mobilization panel

Quick and easy access for service teams


Mattress holders

Bed length extension

Four sockets for mounting accessories

SoftDropPlus™ side rails that lower below the mattress platform

Central wheel lock





# ECONOMY

While designing the NEXO bed we wanted to offer patients comfort and safety and hospital staff reliable equipment, which would be easy to operate and economical. We want you to be able to afford more.

## Legendary Quality and proven effectiveness

The offered NEXO variant is next version of the popular hospital bed, the reference base of which reaches almost one hundred countries around the world. Thousands of units proved the quality and durability, and the feedback from the market allowed for even better development of functionalities. Buying a NEXO means buying a symbol of quality and experience.

## Low service costs

Famed NEXO is a product of the legendary producer of medical devices – Żywiec factory. A Polish product also means Polish service - with decades of experience and a short response time. NEXO was created with the help of our service technicians, so many service activities are carried out in a simple and quick way.

## Increasing biological safety

Fewer places for dirt to collect, a variant for decontamination stations and pSilver™ plastics and varnish technology help the hospital in the work against the risk of infection. This reduces the amount of HAI (Hospital Acquired Infections), which shortens the patient's recovery time and reduces the costs of treating those infections.

## Versatility - one bed for many needs

Full compliance with the EN 60601-2-52 standard, full-length side rails, high stability, controls locks, compatibility with orthopedic accessories and functionalities allowing for use in ICU departments ensure high versatility of the product, depending on the current number of patients in the hospital.

## Version for BDS

With a huge investment in the automation of the cleaning and disinfection process, it is equally important to choose a durable and well-thought-out design that has proven its quality in hundreds of washing cycles. The NEXO bed has been adapted to decontamination stations with temperatures reaching 65-80°C.

## More beds for your facility

Nexo will meet the need to improve the patient's comfort without compromising the hospital's level of safety and will enable you the purchase more beds, translating into increased profitability of the facility.

## Reduce re-preparation time for a new patient

The reduced number of items operators have to pull out for cleaning and the reduced number and size of spaces where contamination can enter means shorter cleaning and disinfection times - staff can get more beds ready to go back to work in shorter period of time.



## 65% OF THE ICU PATIENTS CAN BE MOBILIZED OUT OF THE BED<sup>4</sup>

Early mobilization is one of the most discussed topic in terms of the ICU. Study shows that out-of-bed rehabilitation is possible in over 65% of the intensive care patients.

One of the most important thing to consider is the neurological state of the patient - general calmness and ability to follow simple commands are required. An out-of-bed mobilization is any activity where the patient sits over the edge of the bed (dangling), stands, walks or marches.

The correct level of mobilization should be determined by the patient's condition and all the safety criteria<sup>4</sup>:

- respiratory considerations,
- cardiovascular considerations (especially the presence of medical devices),
- neurological considerations,
- other (after surgery requirements/general condition).



The side rails lower below the mattress and below the mattress platform.

This allows patient to slide off the mattress safely when standing up, without the risk of pinching his body. Additionally, it protects the side rail mechanism against breaking.



# VEM

**Very Early Mobilization - Very early mobilization is an attempt to ensure that the patient returns to full health as soon as possible, shortening the time from surgery to the first rehabilitation treatments.**

Very early mobilization of the patient is an essential element of any modern hospital. Famed NEXO makes it easier for the patient to get out of bed and exercise by using dedicated VEM panels, the highest longitudinal tilt parameter on the market, as well as foot board lock, stable boards and side rails with many places for a firm grip. Famed NEXO can help the patient get up and away from bed, encouraging him to regain shape and health. For the hospital, this means faster discharge of patients who are as close to their full recovery as possible.

#1

**First phase** - the patient can use the panel integrated into the siderail to move the backrest up, higher and higher - as long as it is comfortable. For the ICU patients the staff should perform the action - use a 30° position or above to prevent VAP. Using a lower angle can be one of the factors for VAP<sup>5</sup>.

#2

In the **second phase** - assisted by medical staff, automatic cardiac chair position can be activated - a staff member assists. The patient can be verticalized a little bit to test if he is in a condition that will enable him to leave the bed in the third phase. Reverse Trendelenburg position is set to test the ability of standing up. Patient can use the locked bed end as a support area for his feet. Some core exercises can be performed and a physiotherapist can start a muscle strength rebuild program.

#3

In **third phase**, with the assistance of staff (3B) and setting the bed in the reverse-Trendelenburg position (the best 18.5° angle on the market), the patient can try to get up. If the patient is strong enough to try to stand on his own, he can use the dedicated control panel located on the outside of the side rail (3C). The anatomical hand grip on the side rail is equipped with up / down buttons to raise the patient to a level where he can stand up on his own.









# FUNCTIONALITIES - AVAILABLE OPTIONS



Wired remote control for the patient and staff (variant with 8 buttons in the Comfort and Premium version).



Wired remote control with mechanical lock (key), 10-button. Remote control for version without the central panel.



A panel for the blind and partially sighted (3D-profiled, contrasting markings).



Panel for personnel with quick access to the Trendelenburg tilt.



VEM mobilization buttons for the patient (up and down movement available when holding the side rail and standing up).



Central panel (supervisor) with VEM functions, one-button positions and electric CPR.



USB sockets (2.1 A) for charging mobile devices.



Ambient LED light under the bed.



Sockets for mounting the IV poles in four corners of the bed (two sockets for a hand grip are also available in addition to those 4 sockets).



Bed length extension with lock (not available in the Light version) allows for a comfortable stay also for taller patients.



Bedpan holder



Urine bottle holder



Foot board with lock.



Fifth wheel makes it easier to guide the bed in narrow places.



Integrated storage space for the power cord, e.g. when transporting the bed.



Plastic hooks that can be removed and moved.



The NEXO bed configurations are selected to ensure maximum comfort and functionality in various price ranges.



**NEXO Light** has electric control, making it easier to operate the bed - both by the patient and the staff (using a wired remote control). The bed has comfortable, individually locked, large wheels (150 mm). Staff can use the best reverse-Trendelenburg tilt on the market, easy-to-disinfect mattress platform, pSilver™ anti-bacterial technology and the optional IPX6 version.



**NEXO Comfort** additionally introduces the philosophy of very early patient mobilization (VEM) to the ICU and residential departments, thanks to lockable board and a panel with dedicated buttons (also used for rehabilitation). Additional comfort is also offered by a multitude of dedicated positions on the central panel and a central wheel lock. NEXO Comfort can be adapted to taller patients due to the extension of the mattress platform. The look of the bed is made more attractive by colorful graphics in side rails and boards. As an option, you can also order a foot controller and alarm systems.



**NEXO Premium** introduces to the bed the control in the side rails (from the outside - for the staff and inside - 3D, contrasting buttons for visually impaired) and enables full use of the bed functionality in very early mobilization approach. NEXO Premium is equipped with dedicated patient mobilization panels (VEM).

Function	NEXO Light	NEXO Comfort	NEXO Premium
Panel for patient with aids for visually impaired	-	-	●
VEM rehabilitation panel for third phase of mobilization	-	-	●
Floor light lowering the risk of falls	-	-	●
Lockable board (leg side)	-	●	●
Nurse call integration module (WL-99.65)	-	-	○
8 buttons remote, lock from central panel	-	●	●
Light version remote (bed without ACP), key lock	●	-	-
Central lock	-	●	●
Two additional lock levers on the head side (PL-14)	-	○	○
Individual lock	●	-	-
Length extension	-	●	●
VEM central panel	-	●	●
Graphics in side rails (inserts)	○	●	●
X-ray translucent, easily removable bed pallet	●	●	●
Backrest inclination indicator	●	●	●
Foot control (PL-63.462, up/down movement)	-	○	○
Out-of-bed alarm (WL-99.69)	○	○	○
Unlocked wheels alarm (WL-99.68)	-	○	○
Preparation for decontamination tunnels (IP66)	○	○	○

**KEY:** - NONE | ○ OPTION | ● INCLUDED

### Accessories and options (compatible with every version):

- Height extension adapters (mattres up to 22 cm). [PL-60]
- Fifth wheel. [PK-01]
- Orthopedic frame. [RW-09]
- Oxygen tank holder (vertical). [PL-36.6]
- IV pole. [WK-12, WK-17]
- Bedpan holder. [WL-20.8]
- Urine bottle holder. [WL-19.7]
- Hand grip. [UR-7.2]
- X-ray cassette tunnel in back rest. [PL-18.5]
- Linen and control panel shelf. [PL-67.66]
- Mattress for length extension. [MC-47]
- 12 cm foam mattress. [MC-28.05W1]
- 16 cm foam mattress. [MC-28.05W5]
- Foam anti-decubitus mattress - 12 cm. [MP-01]
- Foam anti-decubitus mattress - variants - on request.
- Active air anti-decubitus mattress - variants - on request.
- USB socket. [WL-99.67]
- Single color inserts. [PL-12]
- Graphic inserts - custom design. [PL-12, extra paid]
- Bed tilt indicator. [PL-24.4]



## Technical data

Length	2270 mm (+ 320 mm bed length extension)
Outer width (with side rails)	990 mm
Internal width (between side rails)	885 mm
Mattress width	850 / 860 mm
Max. mattress thickness	160/220 <sup>6</sup> mm
Height reulation (min - max)	390 - 800 mm
Maximum backrest angle	70°
Maximum thigh segment angle	40°
Trendelenburg tilt	16,5°
Reverse-Trendelenburg tilt	18,5°
Backrest autoregression	120 mm
Autoregression of the thigh section	45 mm
Clearance under the bed (height x length)	150 x 1500 mm
Percentage of chassis cove by the mattress platform	91%
Ingress protection (according to IEC 60529 standard)	IPX4 / IPX6 (option)
Maximum load (SWL)	250 kg
Wheel diameter	150 mm



<sup>1</sup> HOSPITAL-ACQUIRED CONDITIONS, DATA BY HOSPITAL ON NATIONALLY STANDARDIZED METRICS - The Leapfrog Group

<sup>2</sup> Guidelines for prevention of hospital acquired infections, Indian Journal of Critical Care Medicine, 2014 Mar; 18(3): 149-163

<sup>3</sup> WHO guidelines on hand hygiene in health care: A summary. 2014. Mar 10

<sup>4</sup> Early Mobilization of Patients in Intensive Care: Organization, Communication and Safety Factors that Influence Translation into Clinical Practice, Critical Care volume 22, Article number: 77 (2018)

<sup>5</sup> Hunter J.D., Corry P.R.: „Ventilator – associated pneumonia” BJA CEPD 2002

<sup>6</sup> Requires the addition of side rails height extension - PL-60.



Famed Żywiec Sp. z o.o.

Fabryczna 1  
34-300 Żywiec  
Poland

Biuro:

tel.: +48 33 866 63 75  
fax: +48 33 475 58 90

export@famed.com.pl  
www.famed.com.pl



The product meets the requirements of the European Directive MDD 93/42/EEC on medical devices and the Act on Medical Devices.

Edit 184/2020/11/1. Famed Żywiec reserves the right to modify the product and specifications along with technical progress. All illustrations, parameters and photos used in this material are used only for demonstration purposes and may not reflect the final product

